

PRODUKTINFORMATIONEN

Lassen Sie sich 100% kostenfrei
und herstellerunabhängig von uns beraten:



09681 796910



ANFRAGE@4MEDIC.DE



Ihr Kontakt zu uns

09681 796910

info@4medic.de | www.4medic.de
Gerberstraße 11 · 92670 Windischeschenbach

4 { m e d i c TM

Q7

Handle Your Care



esaote
HEALTH WITH CARE

Q7

Handle Your Care

The Q7 wireless ultrasound system takes convenience and portability to the next level.

Its ultra-small size and light weight allow you to carry it right in your pocket, making it ideal for on-the-field use.

This affordable and user-friendly device is the perfect tool to enhance patient outcomes and support your diagnosis.



ULTRA-COMPACT DESIGN

Hyper-light weight (160 g) construction and compact design (150x40x25 mm) for on-the-fly scanning.




6 STREAMLINED BUTTONS

Facilitating control and image optimization for a simplified workflow.



STABLE WIRELESS CONNECTION

Support for Windows®, Android® and IOS® operating systems on PCs and smart devices through secure data flow management.



FAST START

Ultra-fast boot-up time (less than 5 seconds) for a quick response in any situation.



CLEANABILITY & DURABILITY

Fully waterproof for quick and easy cleaning. Built to last in any environment.



EXTENDED BATTERY LIFE

More than 4 hours of high-performance scanning with low power consumption (3W) to increase your productivity.



FLEXIBLE DATA MANAGEMENT

Flexible wireless dataflow with DICOM PACS and Worklist server query. User authentication and auto-anonymized images for safe handling of clinical information.



3 YEARS WARRANTY




Essential protection to keep your investment running safely.



Quality Image at Your Fingertips

Q7 is easily adaptable to many ultrasound examinations and clinical applications thanks to its interchangeable technology, which allows you to change only the transducer head while leaving the base in place.

To improve your workflow and save even more time, find out more about the dual-use heart&ABS probe, which is designed for routine cardiac and abdominal examinations.





160000539 Ver.01



Youkey, 6th F, Block B, No. 1 Building, B1 Area, Hantian Tech City, Nanhai District, Foshan City, Guangdong, China

Imported and distributed by Esaote S.p.A. - sole-shareholder company

Via Enrico Melen 77, 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275,
info@esaote.com - www.esaote.com

Android® is a registered trademark of Google. Windows® is a registered trademark of Microsoft Corporation. IOS® is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license. Technology and features are device/configuration-dependent. Specifications subject to change without notice. Information might refer to products or modalities not yet approved in all countries. Product images are for illustrative purposes only. For further details, please contact your Esaote sales representative.

Please visit us online
for more information

